




- Combining the advantages of durable and disposable
- Flexibility in the use of cages
- Cost effective solutions for contingency planning



HOUSING

TECNIPLAST
innovation through passionSM



HOUSING

Cage Liner Disposable Cage Insert

1/ Introduction: disposable cage liners

Tecniplast has recently developed a series of liners to be used inside existing plastic cage bodies.

Those are made of a very thin layer of Polyethylene stackable one inside the other, 100 per crate.

2/ Combining the advantages of durable and disposable

Tecniplast disposable cage liners are designed to be used in durable plastic cage, combining the solidity and reliability of a durable plastic cage with the advantages of a disposable one.

3/ Flexibility in the use of the cages

The liners have been developed in order to address unmet market needs for disposable cages for special applications such as:


- Irradiation studies
- Imaging
- Toxicology studies
- Quarantine
- Nuclear Medicine studies (isotopes' treated mice)

4/ Cost effective solutions for contingency planning

Disposable liners can be used in specific situations, in contingencies where washing capabilities are limited or non-existent.

Disposable cage liners

1



Flexibility in the use of the cages

2&3



TECHNICAL SHEET

GM500DP

GR900DP

1145T00DP

1284L00DP

1285L00DP

1290D00DP

1291H00DP

Plastic: Polyethylene (A-PET).

Thickness: 0.1 ± 1.6 mm.

FDA approved

TECNIPLAST®

HEADQUARTERS

ITALY | Tel. +39 0332 809711 • www.tecniplast.it • E-mail: tecnicom@tecniplast.it

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. +61 2 8845 6500 • www.tecniplast.it • E-mail: info@tecniplast.com.au

FRANCE | Tél. 04 72 52 94 41 • www.tecniplast.fr • E-mail: info@tecniplast.fr

GERMANY | Tel. 08805 921320 • www.tecniplast.de • E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. 0845 0504556 • www.tecniplast.it • E-mail: info@tecniplastuk.com

JAPAN | Tel. +03 5770 5375 • www.tecniplast.it • E-mail: info@tecniplastjapan.co.jp

USA/CANADA | Toll Free: 877.669.2243 • www.tecniplastusa.com • E-mail: info@tecniplastusa.com

To find your local distributor, please visit www.tecniplast.it



ISO-14001 - Cert. n° 0256A/0



ISO-9001 - Cert. n° 0875

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.
Cage Liner is a Tecniplast trademark.

CAT105/ING - Rev. February 2011