# TEC 2000 Paraffin Dispenser

IVD CE



## The unit

TEC 2000 paraffin dispenser is a unit designed for paraffin of different melting points and purpose.
The upper reservoir can contain up to 5 lt of paraffin.
The particular structure of the reservoir allows homogeneous heat dissipation and prevents thermal shock to the molecular structure of the paraffin, which retains its characteristics over time.
Three different areas with settable temperature: wax, tap and working area.
Easy-to read display allows viewing of the parameters at a glance.
Economic, reliable, ideal for small labs.




## Technical features

- Nominal supply voltage:220 V/50 Hz or 110 V/60 Hz
- Maximum power draw:600 V/A
- Working temperature: from 30° to 90°C increments 1°C
- Dimensions: Width: 31 cm, Depth: 35 cm, Height: 41 cm
- Heat plate: width: 31 cm, depth: 15 cm
- Weight: 13 kg

| Product | Size | Packaging | Code |
|---|---|---|---|
| Paraffin Dispenser | Unit | 1 pc | TEC2000 |