TEC 3000 Histo Trimmer

IVD CE



The unit

TEC-3000 Wax Trimmer is a smart, practical device designed to remove excess of paraffin from the embedding cassettes after inclusion of specimens.

Scraping the embedded block on the "grooved" plane heated surface, excess paraffin will be removed leaving the cassette clean and ready to be set properly on the clamp of the microtome.

The excess of melted paraffin is collected in a practical removable drawer.

No more residual paraffin curls on workbench!

Technical features

- Temperature range: from 30° to 90°C increments 1°C
- Precision in the control of temperature: ± 1°C
- Nominal supply voltage: 220 V ± 10% 50 Hz
- Or 110 V ± 10% 50/60 Hz



Dimensions

- Width: 22 cm
- Depth: 16 cm
- Height: 17,5 cm
- Weight: 2,5 Kg

Product	Size	Packaging	Code
Histo trimmer	Unit	1 pc	TEC3000