Product technical data sheet

C41

Cod. POR0041



Disclaimer:
Images are purely indicative. The features and the accessories of the products may differ according to model and application.

Characteristics

Temperature range	+10/-10 °C
Capacity	41 lt (1,45 cu ft)
External Dimensions (WxDxH)	cm 65x38,5x40 *Internal dimensions (WxDxH) cm 37x28,5x30* *(+ additional compartment 15x28,5x12,5); Weight kg 17*
Structure and Insulation	Atoxic plastic materiai int./est., thickness 40 mm
Door	Nr 1 blind lid
Key Lock	
Internal equipment	Nr 1 open wire basket dim. cm. 28x25x25 (WxDxH) Nr.2 cables: far 12/24V and far 110/240V
Controller	Digital electronic thermostat
Alarms	Acoustic and Visual far temperature deviations
Refrigeration Type	Static
Defrost	Manual

Technical data

Details	• Refrigerant type: R600a CFC free • Mains Voltage 12V/24V; 110V/240V/50Hz-60Hz • Nominal consumption 12 Vcc: 2,58 A • Nominal power consumption: 31 W • Ambient Temperature: +10/+32 °C; humidity 60% • Connectors: 2 Sockets (1 for 12/24V and 1 for 230V)
In accordance with	2014/30/CE, 2014/35/CE and following modifications
Package	Strong Cardboard on wooden pallet *Dim. cm.64x42x56; Kg.23* *• Wooden crate on wooden pallet avai/able: dim. cm.75x50x62; Kg.36*


Accessories

Datalogger


FIOCCHETTI SCIENTIFIC S.R.L.
Headquarter: Via Panagulis, 48 - 42045 Luzzara (RE) Italy
VAT No. IT00934960352
Phone +39 0522 976232 - Fax +39 0522 976028
www.fiocchetti.com - Email: info@fiocchetti.it - Pec: fiocchetti@pec.it

Misure interne Frigoriferi portatili/ Portable refrigerators inner dimensions

C29 misure in mm/ dimensions in mm



C41 misure in mm/ dimensions in mm



C65 – C65T misure in mm/ dimensions in mm

