GENESIS™
CENTRIFUGE BALANCE

MODEL# CB220
CAT# 419-CB220



Now it's easy to keep things in balance.

- Electronic differential digital scale
- Displays weight differential of 2 objects
- Balances in seconds
- Simple and easy to operate
- Protects expensive centrifuges from unbalanced loads

GENESIS BPS™

Technology is in our Blood.℠

1-866-71-BLOOD | www.GenesisBPS.com

GENESIS™
CENTRIFUGE BALANCE

MODEL# CB220
CAT# 419-CB220



SPECIAL FEATURES

- Electronic differential digital scale
- 2 individual weight sensors and control circuitries
- Independently displays weight differential of 2 objects
- Indicator lamp displays which object weighs less and by how many grams
- Balances in seconds
- Simple and easy to operate
- Protects expensive centrifuges from unbalanced loads

ACCESSORIES INCLUDED WITH PURCHASE OF GENESIS™ CENTRIFUGE BALANCE

- A/C Power adapter
- Spare Fuse 250V/1A
- User Manual

SPECIFICATIONS

Dimensions ..... 13"L x 7"W x 5"H
..... 330L x 180W x 130H mm
Weight measurement ..... 0-2500 grams (range)
Maximum weight difference ..... 2500 grams
Weight sensors ..... Load cells
Weight display panel ..... 4-digit LED display

INDICATOR LAMPS

Good (if both trays are balanced)
A (left side of balance tray)
B (right side of balance tray)

ACCURACY

Weight calibration tolerance ..... +/- 1 gram
Decimal truncation error ..... 1 gram
Total accuracy of display ..... +/- 2 grams
Weight accuracy ..... 4 grams

POWER SOURCE

18V DC (A/C power adapter 100-120/220-240V A/C in)
Power consumption ..... 5 watts

WEIGHT

Centrifuge Balance ..... 7.1 lb
Shipping weight ..... 17.1 lb
Package ..... one unit per package

To order, call Genesis BPS at

1-866-712-5663

GENESIS BPS™

Technology is in our Blood.℠

Genesis BPS | 65 Commerce Way, Hackensack, NJ 07601 | TEL: 866-71-BLOOD | FAX: 201-708-1104 | info@GenesisBPS.com | www.GenesisBPS.com