ALCHIMIA

PSS-L

PHOSPHATE BUFFERED SALINE SOLUTION

SOLUTION FOR DONOR GLOBE AND CORNEAL RINSING

GLOBE AND CORNEAL RINSING

	PSS-L
Ingredients	‣ sodium chloride ‣ sodium phosphate ‣ potassium chloride ‣ potassium phosphate ‣ sodium hydroxide
Form	sterile liquid
Volume	30 ml
Physico-chemical parameters	‣ hydrogenionic activity (pH): 7.20 - 7.60 pH units ‣ chloride content: 0.80 - 0.88 %w/w ‣ osmolality: 255 - 345 mOsm/Kg
Sterilization	aseptic filtration and filling according to UNI EN 556-2 and UNI EN ISO 13408-1
Storage temperature	+2°C / +25°C
Expiry date	48 months
Intended use	donor globe and corneal rinsing
Temperature of actual use	room temperature
Class	medical device
GMDN code	P36094
CND code	Q020299
MD registration number	3736/R
Medical device class	IIa
CE mark	CE 0546

PSS-L is a phosphate buffered saline solution for donor globe rinsing, before corneal excision, and for corneal rinsing during trypan blue evaluation.

Ref.: GRS 003-00

Contents:

- 12 30ml-vial pack
- 1 package insert

GS1-128 barcode

ALCHIMIA

Viale Austria, 14
35020 Ponte S.Nicolò (PD)
Italy

VAT IT 00063370282

Phone +39 049 8962074 / 64
Fax +39 049 8962071

www.alchimiasrl.com
info@alchimiasrl.com