

Quality at a reasonable price

PARI COMPACT

- Inhalation device for diseases of the lower airways, such as **COPD, Asthma or Bronchitis**
- **Proven efficacy:** includes PARI LC PLUS® nebuliser with mouthpiece and child's mask
- Can be used in combination with the PARI LC SPRINT® XL*ent* nebuliser to treat diseases of the upper airways such as croup, laryngitis, pharyngitis.

Technical data:

Operating pressure: 1.0 bar
Weight: 1.2 kg
Dimensions (w x h x d) 10 x 14 x 16 [cm]
Warranty: 2 Years

Aerosol Performance with PARI LC PLUS®

Total Output Rate: 418 mg/min
MMD: 3.9 µm (Mass Median Diameter)
Mass percentage below 5 µm: 64 %

Aerosol Performance with PARI XL*ent*

Total Output Rate: 680 mg/min
MMD: 8.5 µm (Mass Median Diameter)
Mass percentage below 5 µm: 68 %

Measured with the Malvern MaserSizer X (according to Fraunhofer-Model) at 21°C and 50% relative humidity.
Nebulised medium: 0.9% NaCl (5 ml). Inspiratory flow 20 l/min

Specialists in effective inhalation PARI

www.pari.de

PARI COMPACT

and its accessories



PARI COMPACT

Inhalation device
230–240 V~/50Hz
With PARI LC PLUS® Nebuliser, Child Mask, Mouthpiece universal and Connection tubing

Item No. 052G1002



Connection Tubing (f/f)

Item No. 041B4592



PARI LC PLUS® Year Pack

consists of LC PLUS® Nebuliser, Connection Tubing and Filter

Item No. 022G8010



PARI Filter Holder with Filter

Item No. 041B6500



PARI Child Mask Soft „Spiggy"

with Elastic Band

Item No. 041G0741



PARI Adult Mask Soft

with Elastic Band

Item No. 041G0470



PARI LC SPRINT® XL*ent*

consists of PARI LC SPRINT® XL*ent* Nebuliser, Connection tubing and Filter

Item No. 023G1801



PARI COMPACT
Item No. 052G1002

PARI
Export Service Center
+49 (0) 8151 279-220

Give us a ring ...

PARI GmbH
Specialists in effective inhalation
Moosstrasse 3 · D-82319 Starnberg
Telefon +49 (0) 8151 279-0
Fax +49 (0) 8151 279-101

You can find PARI on the internet:
E-Mail: info@pari.de · **www.pari.de**

„misprint excepted"

Specialists in effective inhalation

© PARI GmbH · Spezialisten für effektive Inhalation · 052D0017-B-09-2014